

Smart Use Smart Share

RIOSensor





Your Clinic will be closer to SMART Digital Dentistry through the combination of RIOSensor and Windows tablet PC. It will improve efficiency of workflow and consultation.

Smart Use

Now, Share Intraoral X-ray images through any kind of smart phone, tablet, PC and laptop with RAYSCAN web providing necessary functions for most Intuitive user interface.



Smart Share

RIS500

Specifications

Outside dimensions
Image area
Pixel dimensions
Sensor Technology
Scintillator
Number of pixels
Image Resolution
Connection
Cable length

Size 1

39x25mm
30x20mm
20x20um
Advanced CMOS with optical fiber
CsI. or Gd₂O₂S
1.5 megapixels
25 lp/mm
USB 2.0 high speed
2 meters



Size 2

42x30mm
34x26mm
20x20um
Advanced CMOS with optical fiber
CsI. or Gd₂O₂S
2.2 megapixels
25 lp/mm
USB 2.0 high speed
2 meters

